UX10 Exam Chair

The UX10 Exam Chair's breakthrough design maximises everyday patient comfort and caregiver efficiency for examinations and treatment.



Key Benefits

- | Three (3) Sections: Full lie-flat to chair positioning
- | Storage points for control and equipment
- | Built in backrest push handles for ease of movement
- | Singular formed base cover
- | Out of the way control and power cable management
- | Forward tilting and side egress
- | Designed, manufactured and assembled in Australia

Application

A specially designed examination chair and bed in one – for new levels of comfort, safety, ease and flexibility, even in the tightest spaces.

Why Buy?

We designed and developed it with specific input by one of Australia's largest metropolitan hospitals – ensuring every detail has practical day-to-day benefits in today's dynamic healthcare environment. Instantly adjustable from chair to lie-down table, its electronic height variability takes ease and safety to new levels. Patients can sit upright, recline, lie flat or even tilt – with the stability of twin powered columns under the chair. They can move on and off the chair with confidence, with its unique ultra-low height option. The lightweight design and large castors make it a breeze to move around consultation rooms, clinics or hospital departments.

Assembly Guide.

v-13052024



UX10 Exam Chair



Ultra-low height

Make it easier and safer for patients to get on and off.

The very low seat height of this chair with fully retracted sides helps with easy egress and takes away the requirement to use step stools or unnecessary heavy lifting. Sit and position your patient comfortably and then raise them to the desired work height with the touch of a button.



Lies flat

The backrest and leg rest move independently or together from chair to bed.

You have the flexibility of sitting a patient up to almost 90 degrees or lying down completely flat. There is a CPR function to lie the backrest flat quickly in case of emergency.



Head and foot tilt

Up to 14-degrees trendelenburg positioning with tilt towards patients' head or feet.

Transition a patient into the trendelenburg position quickly in an emergency.



Easy transfer

Conveniently positioned handles and large wheels for ease of location.

Lock in the steer castor and use the push handles to move the chair effortlessly to where you need it.



Low - High seated

Great low - high range for all medical applications.



Low - High lie flat

Patient friendly low - high lie flat positioning.

Specifications

Size: 650mm wide x 1350mm overall (1900mm

lying flat).

Comfort: 75mm thick premium, sealed,

removable cushion.

Height Range: 500mm - 880mm with lie-flat and

forward / reverse tilt.

SWL: 220kg.

Castor Size: Large twin wheel 125mm castors

with 4 x brake and steer.

Product Code

UX10-EC

v-13052024

