

# Gel Pad Patches

CleanPatch®-P repair patch for gel pads and positioners and restores damaged surfaces to an intact and hygienic state in seconds.

## Key Benefits

- Constructed of ultra thin, clear and tough material
- Conforms to curved pressure care surfaces
- Glue seals right to the edges preventing fluid ingress and lift
- Cost effectively extend the life of specialist equipment
- Fully cleanable with hospital grade disinfectants
- Clinically validated and TGA registered

## Product Code

- 10 pack - CPGP-31001

## Patch Dimensions

- 45 x 89mm rectangle patch
- There must be a >25mm margin around the damage site to ensure optimal adhesion.

## Why Buy

The gel pad patch has been designed to repair surfaces of gel positioners commonly used in operating room environments. The rectangle shaped patches are very useful for damage such as tears, punctures and abrasion that can be found on these high-risk surfaces. A cost-effective way of quickly ensuring that damage is checked and doesn't remain an infection risk.

[View how to correctly apply the Gel Pad Patch.](#)



v-16042024

"Confidence & Care, every day"

[selectpatientcare.com.au](http://selectpatientcare.com.au)

**select**